ParaPost® Fiber WHITE™

Esthetic Post System

Instructions for Use

CE 0344

 COLTENE

ParaPost®Fiber WHITE™

ParaPost®Fiber WHITE™

EN	Instructions For Use.....	4
DE	Gebrauchsanleitung	8
FR	Mode d'emploi	12
ES	Instrucciones de uso.....	16
IT	Istruzioni per l'uso.....	20
NL	Gebruiksaanwijzing.....	24
SV	Bruksanvisning	28
DA	Brugsanvisning.....	32
FI	Käyttöohjeet	36
NO	Bruksanvisning	40
PT	Instruções de utilização	44
EL	Οδηγίες χρήσης.....	48
LV	Lietošanas instrukcija.....	52
ZH	使用说明	56

One-Office-Visit Technique

Intended Use: Fiber posts are intended to be cemented into the root canal of a tooth to stabilize and support a restoration.

 <p>CAUTION!</p>	<p>Federal law (USA) restricts this device to sale by or on the order of a healthcare professional.</p>
 <p>WARNING!</p>	<p>Failure to properly follow the directions in this IFU may result in patient injury or compromise the integrity of the restoration.</p>
 <p>NON STERILE</p>	<p>Posts are delivered non-sterile and are to be sterilized prior to use.</p>
 <p>WARNING!</p>	<p>Material is composed of a metacrylate-based glass fiber reinforced resin. Device should not be used in patients with known material sensitivity.</p>

Note: Posts are single use devices and should never be salvaged for re-use. Structural integrity of the post could be compromised including microscopic damage which could lead to post fracture.

For illustrations see page 59.

We strongly recommend the use of Hygenic® Rubber Dam with this procedure.

1. After endodontic therapy has been completed, prepare the tooth for restoration as if pulp were intact. The preparation should include at least 1.5mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect (**Fig. 01**).^{*} If that amount of tooth structure is not available, it should be created by periodontal crown lengthening and/or orthodontic extrusion.^{*}
2. Use a radiograph to determine the appropriate diameter and depth of post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1mm minimum)^{*} to prevent perforation or weakening of root walls, and at least 4 - 5mm of gutta percha must be left intact to protect the apical seal.
3. Remove gutta percha to preplanned depth with a Gates-Glidden drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a silver cone, remove it and reseat with Hygenic gutta percha before preparing the post space.
4. To begin paralleling the post space, select the Para-Post Drill diameter which corresponds to the last Gates-Glidden drill used in Step 3 (**Fig. 01**). (See "Approximate Diameter Comparisons" chart, **Fig. A**). To help maintain the preplanned depth of the preparation, ParaPost Drills have been pre-marked to identify depths of 7, 9 and 11mm from the apical end of the drill.

ParaPost Drills can be used either manually with the Universal Hand Driver or with a slow-speed contra-angle (750 - 1,000 RPM). When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the risk of the drill seizing in the post space. Remove all weakened or unsupported tooth structure. Remove any tooth debris from canal by irrigating post space with water spray.

5. Sequentially step up to the next larger ParaPost Drill until the preplanned diameter and depth are achieved (**Fig. 02**).
6. Use a cylindrical diamond or carbide bur to prepare an anti-rotational box (**Fig. 03**).
7. Select the ParaPost Fiber White Post that corresponds to the last drill used to prepare the post space. Remove colored band from post head. Insert post into post space (**Fig. 04**). Check carefully for occlusal clearance. Remove post; shorten as necessary from either the apical or occlusal end (as clinical judgement dictates) with a Carborundum separating disk (**Fig. 05**). (Wet the post with water before cutting and rotate the post when cutting to ensure an even cut of each fiber.)
8. After try-in and before cementing, gently clean post with an alcohol wipe.
9. (Optional) Brush silane onto the post and allow to dry.

10. Condition the canal and/or apply primer/adhesive as per cement manufacturer's instructions (we recommend **ParaPost Cement**).
11. Use the cement components according to manufacturer's instructions and apply to post surface. Use a Lentulo Spiral to coat the walls of the post space with cement, if setting time of cement will allow (**Fig. 06**).
12. Seat post immediately, slowly inserting the post to full depth, allowing excess cement to vent (**Fig. 07**). Apply pressure for about 60 seconds.
13. In the event light curing is required, wipe away any excess cement before activating the light and follow manufacturer's instruction.
14. Complete the core and fabricate the final restoration (**Fig. 08**).

Cleaning and Sterilization

Posts: Initiate cleaning by soaking the post(s) for 5 minutes in a solution consisting of an enzymatic cleaner, such as Terg-a-Zyme and de-ionized water at approximately 20°C (68°F). Follow the agent manufacturer's instructions.

Clean the post(s) by brushing away from the post while keeping it immersed in an enzymatic cleaning/disinfecting agent validated for use on dental instruments. Special care should be taken

to clean crevices and other hard to reach areas. Thoroughly rinse with tap water. Dry using a paper towel.

Additional cleaning can be performed in an ultrasonic cleaner for 5 minutes. Remove and thoroughly rinse with tap water; dry with a paper towel.

Place post(s) into a small, clean covered glass container.

Add a ready-to-use liquid sterilant solution containing 7.50% hydrogen peroxide and 0.85% phosphoric acid, (Sporox II**), covering the post(s). Let stand at room temperature for a minimum of 6 hours but not to exceed 24 hours. Pour off solution into sink and add sterile water to glass container covering the post(s). Gently agitate rinse solution for two minutes, then discard solution into sink and dump post(s) onto sterile sponge in clean container. Affix cover and store until ready for use.

Follow solution manufacturer's instruction.

Never re-use solution - never add fresh solution to used solution.

Always discard solution after each use!

Autoclave Sterilization Option: ParaPost Fiber White may be dynamic-air-removal steam sterilized (4 minutes at 132°C followed by drying for 20 minutes). Follow manufacturer's recommendations. Posts must be sealed in a sterilization pouch

before use in autoclave. Alternatively, sterilize in a pre-vacuum chamber at 134°C for 3 minutes or 134°C for 6 minutes.

Drills: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills, with a sterilization indicator, into an autoclave bag and seal. Sterilize in gravity at 132°C for 12 minutes. Sterilizing bags should not touch the walls of the autoclave. Alternatively, sterilize in a pre-vacuum chamber at 132°C for 4 minutes, 134°C for 3 minutes or 134°C for 6 minutes. Once the sterilization monitor indicates that sterilization was achieved, place autoclave bag, with the drills inside, in storage until ready for use. After use, to decontaminate, soak in a 70% isopropyl alcohol solution for 10 minutes.

Universal Hand Driver: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by autoclave (see "Drills" for autoclave instructions).

Warning Notice: DO NOT BEND POSTS

Post Removal Technique

1. Take a radiograph of the tooth that requires post removal.
2. Carefully estimate length and diameter of ParaPost Fiber White post to be removed.
3. Remove any restorative or core material which may be in the way to gain direct straight-line access to the post.

4. Indent the center of the post using a # 1/2 round bur at high speed.
 5. Carefully drill an initial channel into the center of the ParaPost Fiber White post using a Kodex twist drill (cat. # K95 or K97) at slow speed.
 6. Carefully drill to full depth through the center of the post using an end cutting Tenax Starter Drill (cat # TEDC1, Ø .036"/0.90 mm) or Tenax Drill (cat #TEA11, Ø .043"/1.1mm). Check radiographically at stages.
 7. Gradually work up to original diameter using standard ParaPost drills.
 8. Re-treat endodontically, if necessary.
 9. Re-restore with new post/ core and restoration. Always verify steps via radiograph.
- * References and/or literature available upon request.
- ** Sporox II is a registered trademark of Reckitt & Coleman, Inc.

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Coltène/Whaledent Inc.
Coltène/Whaledent AG
Coltène/Whaledent GmbH + Co. KG