ParaPost XP System

No-Ox Precious Alloy / G-Post Precious Alloy
Cast-To Gold Posts

Intended use:

Posts are intended to be cemented into the root canal of a tooth to stabilize and support a restoration.

Casting Technique – Direct and Indirect

1. After endodontic therapy has been completed, prepare tooth for restoration as if pulp were intact. The preparation should include at least 1.5 mm of sound tooth structure all around the circumference of the preparation, apical to the core, for desired ferrule effect. If that amount of tooth structure is not available, it should be created by periodontal crown lengthening and/or orthodontic extrusion.

2. Use a radiograph to determine the appropriate diameter and depth of post space preparation. Keep in mind that sufficient root wall thickness must be maintained (1 mm minimum) to prevent perforation or weakening of root walls, and at least 4-5 mm gutta-percha must be left intact to protect the apical seal.

3. Remove gutta-percha to preplanned depth with a Gates-Glidden Drill, Peeso reamer and/or hot instrument. Radiographic verification is recommended. If the canal has been obturated with a silver cone, remove it and reseat with Gutta-Percha before preparing the post space.

4. To begin paralleling the post space, select the ParaPost Drill which corresponds to the last Gates-Glidden drill used in Step 3 (Fig. A). To help maintain the preplanned depth of the post space, ParaPost Drills have been premarked to identify depths of 7, 9 and 11 mm from the apical end of the drill. ParaPost Drills can be used either manually with the Universal Hand Driver or with a slow-speed contra-angle (750-1000 rpm). When a contra-angle is used, the drill must be kept in continuous clockwise rotation until it has been completely removed from the tooth. This will minimize the possibility of the drill seizing in the post space. Remove all weakened or unsupported tooth structure. Remove any tooth debris from canal by irrigating post space with water spray.

5. Sequentially step up to the next larger ParaPost Drill until the previously determined diameter and depth are achieved. 6. Use a cylindrical diamond or carbide bur to prepare an anti-rotational box without undercuts.

Direct Only

7. Place into the prepared post space either a ParaPost Burnout Post or a ParaPost No-Ox Alloy / G-Post Alloy Cast-To Gold Post, that corresponds to the largest ParaPost Drill used to pre-

pare the post space. If necessary, shorten the Burnout Post or No-Ox / G-Post Alloy Cast-To Gold Post from the apical end.

8. Build up the core pattern with a pattern resin.

9. Remove the post/core pattern from the tooth.

10. Insert the corresponding ParaPost Temporary Post into the prepared post space. If necessary, shorten from the apical end. Fabricate a provisional crown with self-curing resin. To preserve the integrity of the post space, apply temporary cement ONLY to the margins of the crown; DO NOT PLACE any cement into the post space or on the post.

11. Sprue, invest and cast. A. When using a ParaPost Burnout Post, cast with any dental alloy. B. When casting to a ParaPost No-Ox Alloy / G-Post Alloy Cast-To Gold Post, use any dental alloy with a casting temperature no higher than 2500 °F (1370 °C).

CAUTION: DO NOT USE gypsum-bonded investment materials or those containing chloride, as they will make the prefabricated post brittle. Suggested investment materials include: Ceramigold, WM80, Ceramigold 2, Cristobalite Inlay, Hi-Temp, Cristobalite Model, Hi-Temp 2, Hydrovest, Chrome Investment X20.

12. Finish and cement the cast post and core.

13. To prepare for cementation, clean and rinse the canal and dry with paper points. Use the dental cement of your choice according to the manufacturer's instructions.

14. Use a Lentulo spiral to coat the walls of the post space with cement.

15. Coat the post shaft with cement and slowly insert the post to full depth, allowing excess cement to vent.

16. In order to obtain maximum retention to crown, DO NOT POLISH the core. Fabricate the final restoration. After the cement has set, remove any excess before completing the core and fabricating the final restoration.

Indirect Only

7. Place into the prepared post space the ParaPost Impression Post corresponding to the last drill used. If necessary, adjust the length of the Impression Post so the post does not contact the impression tray.

8. Take an impression with an elastomeric impression material (e.g. Affinis). Make sure the anti-rotational box is filled with impression material. NOTE: To ensure retention of the post in the impression material when using Hydrocolloid, form a nail head on the top of the Impression Post by pressing on the head with a hot instrument. DO NOT ALTER the plastic Impression Post when using silicone or rubber-based impression material.

9. Insert the corresponding ParaPost Temporary Post into the prepared post space. If necessary, shorten from the apical end. Fabricate a provisional crown. To preserve the integrity of the post space, apply temporary cement ONLY to the margins of the crown; DO NOT PLACE any cement into the post space or on the post.

10. Pour a model and fabricate a precision removable die using the Pindex® or Accutrac® System.

11. Place the corresponding ParaPost Burnout Post or ParaPost Precious No-Ox Alloy / G-Post Alloy Cast-To Gold Post into the post space of the model. If the No-Ox / G-Alloy Cast-To Gold Post does not readily fit the post space in the die, carefully remove the post space with the corresponding ParaPost Drill in the Universal Hand Driver. If necessary, shorten the Burnout Post or No-Ox Alloy / G-Post Alloy Cast-To Gold Post from the coronal end.

12. Complete the core pattern.

13. Carefully remove the post/core pattern.

14. Sprue, invest and cast. A. When using a plastic ParaPost Burnout Post, cast with any alloy. B. When casting to a ParaPost No-Ox Alloy / G-Post Alloy Cast-To Gold Post, use any dental alloy with a casting temperature no higher than 2500 °F (1370 °C). CAUTION: DO NOT USE gypsum-bonded investment materials or those containing chloride, as they will make the prefabricated post brittle. Suggested investment materials include: Ceramigold, WM80, Ceramigold 2, Cristobalite Inlay, Hi-Temp, Cristobalite Model, Hi-Temp 2, Hydrovest, Chrome Investment X-20.

15. Finish and cement the cast post core.

16. To prepare for cementation, clean and rinse the canal and dry with Paper Points. Use the dental cement of your choice according to the manufacturer's instructions.

17. Use a Lentulo spiral to coat the walls of the post space with cement.

18. Coat the post shaft with cement and slowly insert the post to full depth, allowing excess cement to vent.

19. In order to obtain maximum retention to crown, DO NOT POLISH the casting. After the cement has set, remove any excess before completing the core and fabricate the final restoration.

Cleaning and Sterilization

Drills: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush; thoroughly rinse and dry. Place drills with a sterilization indicator, into an autoclave bag and seal. Sterilize in autoclave at 132 °C for 12 min. Sterilizing bags should not touch the walls of the autoclave. Once the sterilization monitor indicates that sterilization was achieved, place autoclave bag, with the drills inside, in storage until ready for use. To decontaminate, soak in a 70 % isopropyl alcohol solution for 10 min.

Posts: Prior to each use, place posts, with a sterilization indicator, into an autoclave bag and seal. Sterilize in an autoclave at 130-135 degrees C for 12-15 min, following standard procedures.

Universal Hand Driver: Prior to each use, remove debris ultrasonically or by hand scrubbing with a brush. Sterilize by autoclave or chemical vapour sterilization.

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Coltene/Wahledent AG
Coltene/Wahledent Inc.
Coltene/Wahledent GmbH + Co. KG