


DenFil Etchant

Semi-gel Type Etching Agent with Phosphoric Acid



FEATURES AND BENEFITS

- DenFil Etchant is a semi-gel type etching agent containing 37% phosphoric acid thickened with natural polymeric materials.

INDICATIONS

- Acid etching of enamel and dentin.

PACKAGE

		ALL PRODUCTS	PRODUCT CODE	KFDA Reg. CODE
DenFil Etchant	Bulk Kit	Etchant 5ml x 5Syringes, Disposable tips (Blue)	ET125-B000	DenFil Etchant-37
		Etchant 3ml x 5Syringes, Disposable tips (Blue)	ET115-E000	
	Eco Kit	Etchant 3ml x 12Syringes, Disposable tips (Blue)	ET136-K000	
		Etchant 30ml x 1Syringe, Disposable tips (Blue)	ET130-B000	